




RAMA  
UNIVERSITY

[www.ramauniversity.ac.in](http://www.ramauniversity.ac.in)

FACULTY OF ENGINEERING & TECHNOLOGY  
DEPARTMENT OF BIOTECHNOLOGY

- 3. This mixture is kept in a PCR equipment and is subjected to repeated cycles of DNA denaturation-renaturation-DNA replication.
- 4. During this process, the decaoligonucleotide will pair with the homologous sequence present at different locations in the DNA.
- 5. DNA replication extend the decaoligonucleotide and copy the sequence continuous with the sequence with which the selected oligonucleotide has paired.



- 6. The repeated cycles of denaturation-renaturation-DNA replication will amplify this sequence of DNA.
- 7. Amplification will take place only of those regions of the genome that have the sequence complementary to the decaoligonucleotide at their both ends.



- 8. After several cycles of amplification the DNA is subjected to gel electrophoresis.

Initial denaturation at 94°C for 10 min.




Denaturation at 94°C for 1 min. Annealing at 37-45°C for 1 min. Extension at 72°C for 1 min.

Final extension at 72°C for 10 min.

Cooling at 4°C.

- 9. The amplified DNA will form a distinct band. it is detected by ethidium bromide staining and visible fluorescence's under U.V. light





**RAPD**

Temperature



Time

RAMA  
UNIVERSITY



KANPUR

3' 5'

Heat

**RAPD**

5' 3'





**RAPD**






Temperature



Time

3' 5'




**RAPD**

5'




5' 3'

Temperature



Time

3' 5'



**RAPD**

